Restricted Product – Prescribe 2 week supply / Reassess

Topical Negative Pressure

To access a device please contact the Tissue Viability Team

♠¥V.A.C ® - KCI Medical

Granufoam Small 10x7.5x3.3cm, Medium 18x12.5x3.3cm, Large 26x15x3.3cm KCI NPWT Gauze Dressing Kit (One Size)

Activac Canister

V↑Renasys ® (Smith & Nephew)

Renasys F (Foam) with Softport Small, Medium, Large

Renasys G (Gauze) with Softport Small, Medium, Large

Renasys Go Canister

Wound Drainage Bags

◆↑Eakin Wound Pouches ® Fold & Tuck Small (wounds 45x30mm) Medium (wounds 110x75mm)Large (wounds 175x110mm) Extra-large (horizontal wounds 245x160mm) Bung Closure Small (wounds 45x30mm), Small Plus (wounds 86x60mm) Medium (wounds 110x75mm) Large (wounds 175x110mm) Extra-large (horizontal wounds 245x160mm) Extra-large (vertical wounds 160x245mm) Extra-large (vertical wounds 290x130mm)

Antimicrobials Prescribing notes – Maximum 2 Weeks Use – Only to be used to treat local signs and symptoms of wound infection. Silvers

- ↑Actisorb®Silver 220 6.5x9.5cm,10.5x10.5cm,10.5x19.5cm
- ↑Tegaderm Alginate Ag 5x5cm,10x10cm,3x30cm rope
- ↑Acticoat ® Flex 3 (3 day)- 5x5cm,10x10cm,10x20cm,20x40cm
- ↑Acticoat ® Flex 7 (7 day) 5x5cm,10x12.5cm,15x15cm
- ↑Flamazine ® 20g, 50g

Honey

- ↑Activon Gel 25g (Tube) Tulle 5x5cm, 10x10cm
- ↑Algivon (Honey & Alginate) 5x5cm, 10x10cm

Glucose Based dressings

- ↑Flaminal Forte (Gel & Alginate) 15g, 50g
- ↑Flaminal Hydro (Gel) 15g, 50g

lodine

- ↑Inadine®- 5x5cm,9.5x9.5cm
- ↑lodosorb® 3g Powder, 10g ointment / lodoflex® 5g,10g,17g Paste

PHMB

↑AMD foam - 5x5cm,10x10cm,15x15cm, 10x20cm, 20x20cm, 8.8x7.5cm Fenestrated

Under Direct Supervision of Tissue Viability Team Only

Vibropulse – Cycloidal therapy

PICO - Disposable NPWT

Larvae Therapy – Freerange Kits and Biobags

Juxta Cures - Compression System

Under Direct Supervision of Podiatry Only

Melonin – 5x5cm, 10x10cm, 20x10cm

Hypafix - 2.5cmx10m,5cmx10m,10cmx10m,15cmx10m

Povitulle - 5x5cm, 9.5x9.5cm



NHS Foundation Trust

Wound Management Formulary - Quick Reference Guide

September 2014 - August 2016

District Nursing, Community Podiatry, Children's Nursing Service, Practice Nurses, and Tissue Viability Team

For additional information please contact

Tissue Viability Team

0191 283 1432

Key / Symbols

♠Primary Dressing Only

◆ Secondary Dressing Only

→ Compression

←Retention / Fixation

↑ Combined Primary / Secondary Dressing

Dressing Packs

Softdrape (order via NHS supply Chain) Small EJA045, Medium EJA046, Large EJA047

Pressure Reduction (Prevention / Cat 1 only)

Kerrapro Pressure reducing pads

Sheet 10x10x0.3cm, 10x10x1.2cm, Strip 30x5x0.3cm Heel One Size

Wound Cleansing

Irripod - 20ml (x25)

Prontosan solution (surfactant) - 350ml (single)

apes

←Cliinpore- 1.25cmx5m,2.5cmx5m,5cmx5m,2.5cmx10m

←Opsite Flexifix - 5cmx1mroll, 10cmx1m roll

←Mefix - 2.5x5m,5cmx5m,10cmx5m, 15cx5m

Debridement Pad

Debrisoft - 10x10cm

Hydrogel

♠Activheal hydrogel - 8g

↑Intrasite conformable -10x10cm, 10x20cm

Skin Barrier / Protection

Medi Derma S Film - 75ml aerosol

Medi Derma S Cream - 90g tube

Film Dressing

↑ Tegaderm film - 6x7cm,12x12cm,15x20cm

Film and Pad Dressing

↑ Tegaderm + Pad- 5x7cm, 9x10cm,9x15cm,9x20cm,9x25cm,9x35cm

Fabric and Pad

↑ Softpore - 6x7cm, 10x10cm, 10x15cm, 10x20cm, 10x25cm, 10x30cm, 10x35cm

Waterproof Fabric & Pad

↑ Mepore Ultra 7x8cm,10x11cm,11x15cm,9x20cm,9x25cm,9x30cm

Low Adherent Wound Contact

↑N/A Ultra - 9.5x9.5cm, 19x19.5cm

↑ Atrauman - 5x5cm,7.5x10cm,10x20cm

↑Telfa Clear- 7.5x7.5cm,7.5x10cm,10x20cm

Hydrocolloids

↑ Hydrocoll Thin Film - 7.5x7.5cm, 10x10cm, 15x15cm

↑ Duoderm Signal - 10x10cm,14x14cm,11x19cm Oval

Hydrofibres

♠Activheal AquaFiber - 2x42cm rope

Alginates

↑Activheal Alginate - 2cmx30cm rope

Super Absorbents

↑Zetuvit Plus - 10x10cm,10x20cm,15x20cm,20x25cm,20x40cm

Dressing Pad

↑Zetuvit E - 10x10cm,10x20cm,20x20cm,20x40cm

Non Adhesive foams

- ↑ Activheal Foam Low to moderate exudate 5x5cm, 10x10cm,10x20cm,20x20cm 18x12cm Heel
- ↑ Allevyn Moderate to high exudate 5x5cm, 10x10cm,10x20cm,20x20cm
- ↑↓Allevyn Gentle Moderate to high exudate 5x5cm,10x10cm,10x20cm,15x15cm,20x20cm
- ↑ Aquacel foam Combination Hydrofiber and Foam high exudate only 5x5cm,10x10cm,15x15cm, 20x20cm,15x20cm

Adhesive Foams

- ↑ VAllevyn Gentle border lite Low to moderate exudate 5x5cm,7.5x7.5cm,10x10cm,5.5x12cm,8x15cm, 10x20cm, 15x15cm, Oval 8.6x7.7cm, 15.2x13.1cm
- ↑ VAllevyn Adhesive Moderate to high exudate

7.5x7.5cm,10x10cm,12.5x12.5cm,12.5x22.5cm,17.5x17.5cm,22.5x22.5cm

- ↑ Vallevyn Gentle Border *Moderate to high exudate* 7.5x7.5cm,10x10cm, 12.5.x12.5cm, 15x15cm, 17.5x17.5cm,10x20cm, Sacrum 16.8x17.1cm, multisite 17.1x17.9cm, Heel 23x23cm
- ↑ ♦ Allevyn Life high exudate only- to be used when changing foam dressings more than x3 per week 10.3x10.3cm, 12.9x12.9cm, 15.4x15.4, 21x21cm
- ↑ Aquacel Foam Combination Hydrofiber and Foam high exudate only

8x8cm,10x10cm,12.5x12.5,17.5x17.5,21x21cm,25x30cm

Retention Bandages

- ←Hospilite® 10cmx4.5m, 15cmx4.5m
- ← Hospicrepe® 239 10cmx4.5m, 15cmx4.5m
- ←Profore® #2 10cmx4.5m

Padding

- ←Profore® #1 10cmx3.5m
- ←Flexi-ban- padding ®- 10cmx3.5m

Tubular bandages

- ←Comfifast® Sizes Red 3.5cmx1m,Green 5cmx5m,Blue 7.5cmx5m,Yellow10.75cmx5m
- ←Comfigauz® Not available on FP10 Sizes 00 Toes, 01 fingers and toes, 56 adult limbs, 78 large adult limbs

Compression Bandages

- → Profore® Kits Elastic Compression Profore lite, <18cm, 18-25cm,25-30cm,>30cm
- → Coflex TLC Malodour Control 2 layer kit Lite, Full Compression, Full Compression XL Kit
- → Coflex UBZ 2 Layer Kit (1st Layer is Zinc impregnated)- Coflex UBZ lite kit, Coflex UBZ Kit
- → Actico® Short Stretch 4cmx6m,6cmx6m,8cmx6m,10cmx6m,12cmx6m
- → Coban® 2 Attenuated Stretch Compression Coban 2 lite Kit, Coban 2 Kit
- → Coban® 2 Components- Comfort Layer 5cmx1.2m,10cmx3.5m, 15cmx3.5m-Compression Layer 5cmx2.7m,10cmx4.5,15cmx4.5m
- →K TWO® Latex Free Kits- 18-25cm, 25-32cm, Reduced18-25cm, 25-32cm

Paste Bandages

Viscopaste ®PB7 - 7.5cmx6m

Compression Hosiery

- Class 1 Superficial or early varices and prevention of deep vein thrombosis while travelling.
- Class 2 Medium varices / Treatment and prevention of venous leg ulcers and associated conditions / Mild oedema

Class 3 - Gross varices / Gross oedema / Treatment and prevention of venous leg ulcers and associated conditions /Post-thrombotic venous insufficiency

Irrespective of system ensure you measure limb using manufacturer guidelines. Please complete manufacturer order form to record the measurements and define class required.

This form must be submitted with prescription to the pharmacy or with the dressing request form to the GP

Mediven Elegance® Standard and Petite Lengths

Mediven For Men ® Standard and Petite Closed Toe

Mediven Plus® Standard and Petite Lengths and Standard and Extra Wide Calf

Sizes I, II, III, IV, V,VI, VII

RAL Classification

Class 1 = 18-21mmHq. Class 2 = 23-32mmHq. Class 3 = 34-46mmHq

ActiLymph® Standard and Petite Lengths

Sizes S, M, L, XL, XXL (XXL Class 1 only)

European Classification

Class 1 = 18-21mmHg, Class 2 = 23-32mmHg, Class 3 = 34-46mmHg

Activa Hosiery®

Unisex ribbed and Patterned Sock (Class 1 and 2 only)

Liners Kit (10mmHg)

Hosiery Kit (40mmHg)

Sizes S, M, L, XL, XXL

British Standard Classification

Liners = 10mmHq, Class 1 = 14-17mmHq, Class 2 = 18-24mmHq, Class 3 = 25-35mmHq